

N95 FOLDABLE MASK

MODEL: ZYE-02



*Images for marketing purposes only



Accreditations:

42 CFR 84 NIOSH N95 

ISO 9001:2015

EN ISO 13485:2016

Approved by Centers for Disease Control and Prevention (CDC).

CDC Approval Link: <https://www.cdc.gov/NIOSH-CEL/ApprovalDetails?schedule=84A&approvalNum=9247>

N95 Particulate Respirator Mask:

Vertical Flat Fold Design

Four-point elastic head straps for a superior mask fit.

4-ply structure:

1st Layer 38g/m² PP Spun-Bond Outer Layer;

2nd Layer: 20g/m² N95 Grade Melt-Blown High Filtration Layer;

3rd Layer: 20g/m² N95 Grade Melt-Blown High Filtration layer;

4th Layer: 250g/m² Needle-Punched Non-Woven Inner Layer

Disposable, non-sterile & foldable.

Latex-Free Knitted Textile Elastic

High Protection, Comfort and Fit

Skin Friendly and does not deform in High Temperatures

N95 - Particulate Filter (95% filter efficiency level). Effective against particulate aerosols free of oil; time use restrictions may apply.


Requires fit testing and must be adjusted to face to provide the intended effectiveness.

Carton Quantity: 400 masks (packs of 20)

Carton Dimensions: 78 x 65 x 50.5 cm

G.W.: 5 kg

Box of 20



NIOSH N95 APPROVED
ZYE-02
N95 PARTICULATE RESPIRATOR



135X120X130MM

Xinze Zhongyi Safety & Protection Products Co., Ltd.
No. 121 Yafan Road, Xuzhou City,
Hubei Province, China 433000



THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC *	Protection [†]	Respirator	Cautions and Limitations *
844-9247	N95	ZYE-02 X	ABCJMNDP

1. Protection

N95-Particulate Filter (95% filter efficiency level)
Effective against particulate aerosols free of oil.
Time use restrictions may apply.

2. Cautions and Limitations

A. Not for use in atmospheres containing less than 19.5% oxygen.
B. Not for use in atmospheres immediately dangerous to life or health.
C. Do not exceed maximum use concentrations established by regulatory standards.
D. Failure to properly use and maintain this product could result in injury or death.
E. All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
F. Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
G. Refer to User Instructions and/or maintenance manuals for information on use and maintenance of these respirators.
H. NIOSH does not evaluate respirators for use as surgical masks.

WARNING:

NIOSH APPROVED: N95
95% efficient against particulate aerosols free of oil. This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting disease or infection. Misuse may result in sickness or death.



Fitting Instructions:
to be followed each time respirator is worn.



1. Preshape top and bottom straps before placing respirator on the face.
2. Cup the respirator in your hand, with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.
3. Position the respirator under your chin with the nosepiece up. Pull the top strap over your head resting it high at the top back of your head. Pull the bottom strap over your head and position it around the neck below the ears.
4. Place your fingertips from both hands at the top of the metal nosepiece. Using two hands mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece. Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.
5. Perform a User Seal Check prior to each wearing. To check fit, place both hands completely over the respirator and exhale. Be careful not to disturb the position of the respirator. If air leaks around nose, adjust the nosepiece as described in step four. If air leaks at the respirator edges, work the straps back along the sides of your head. If you CANNOT achieve a proper fit, DO NOT enter the contaminated area. See your supervisor.

Removal Instructions:

See step 3 of Fitting Instructions and cup respirator in hand to maintain position on face. Pull bottom strap overhead. Still holding respirator in position, pull top strap over head and remove respirator.
This respirator contains no components made from natural rubber latex.


NIOSH Approved: N95

At least 95% filtration efficiency against solid and liquid aerosols that do not contain oil.
PROTECTION:

N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply.



NIOSH N95 APPROVED
ZYE-02
N95 PARTICULATE RESPIRATOR



Box of 20

N95 Particulate Respirator



WARNING
This respirator helps protect against certain particles. **Misuse may result in sickness or death.** For proper use, see supervisor.

IMPORTANT:

Before use, the wearer must read and understand these User Instructions. Keep insert for reference.

Use For:

Solids such as those from processing minerals, coal, iron ore, flour, and certain other substances. Liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapors.

Do Not Use For:

Fine spray, oil aerosols, gases, vapors, asbestos or sandblasting. This respirator does not supply oxygen.

Use Instructions:

1. Failure to follow all instructions and limitations on the use of this respirator and/or failure to wear this respirator during all times of exposure can reduce respirator effectiveness and may result in sickness or death.
2. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the requirements of OSHA 29 CFR 1910.134 such as training and fit testing and applicable OSHA substance specific standards. In Canada, CSA standard Z94.4-93 requirements must be met.
3. The particles which can be dangerous to your health include those so small that you cannot see them.
4. Leave the contaminated area immediately and contact supervisor if dizziness, irritation, or other distress occurs.
5. Store the respirator away from contaminated areas when not in use.
6. Dispose of used product in accordance with applicable regulations.

Use Limitations:

1. This respirator does not supply oxygen. Do not use in atmospheres containing less than 19.5% oxygen.
2. Do not use when concentrations of contaminants are immediately dangerous to life and health are unknown or when concentrations exceed 10 times the permissible exposure limit (PEL) or according to specific OSHA standards or applicable government regulations, whichever is lower.
3. Do not alter above or misuse this respirator.
4. Do not use with beards or other facial hair or other conditions that prevent a good seal between the face and the edge of the respirator.

Time Use Limitations:

If respirator becomes damaged, soiled, or breathing becomes difficult, leave the contaminated area and replace the respirator.



Manufacture Date
Expiry Date
Lot No.





ZHONG YI MASK

NIOSH N95 APPROVED

ZYE-02

N95 PARTICULATE RESPIRATOR



FOR SINGLE USE ONLY

Box of 20

N95 Particulate Respirator



This respirator helps protect against certain particles. Wearing may result in sickness or death. For proper use, see instructions.

IMPORTANT:

Before use, the wearer must read and understand these Use Instructions. Do not use if instructions are not followed.

Use For:

Filtering out at least 95% of airborne particles, such as dust, dirt, and pollen, when used as directed. It is not designed to protect against chemical vapors, gases, or odors.

Do Not Use For:

Protecting against toxic gases, vapors, aerosols, or irritants. This respirator does not protect against infectious agents.

Use Instructions:

1. Perform a series of fit checks and inspections on the use of this respirator before use to ensure it is properly fitted and that the respirator is in good working order.
2. Before use, inspect the respirator for any damage or defects. If the respirator is damaged, do not use it.
3. The person wearing the respirator must be properly fitted and must perform a series of fit checks to ensure the respirator is properly fitted.
4. Label the respirator with the user's name and the date of use.
5. Store the respirator in a clean, dry, and well-ventilated area.
6. Dispose of this product in accordance with applicable regulations.

Use Limitations:

1. This respirator does not protect against:
2. Ozone, oil-based mists, aerosols, or vapors.
3. The use of this respirator is limited to the use of the respirator as described in these instructions.
4. The use of this respirator is limited to the use of the respirator as described in these instructions.
5. The use of this respirator is limited to the use of the respirator as described in these instructions.

Time Use Limitations:

The use of this respirator is limited to the use of the respirator as described in these instructions.

Fitting Instructions:

The use of this respirator is limited to the use of the respirator as described in these instructions.





Zhongji Safety & Protection Products Co., Ltd.
No. 33, Shengyuan Road, Xiangtan City, Hubei Province, China 43000
Tel: 86-731-527025 Fax: 86-731-527045
www.zhongjinisk.com



THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC	Protection ¹	Respirator	Caution and Limitations ²
		ZY6-02	
BAA-6047	N95	X	ABCJMNDP

1 Protection

N95-Particulate Filter (95% filter efficiency level)
Effective against particulate aerosols free of oil;
Time use restrictions may apply

2 Caution and Limitations

- A. Not for use in atmospheres containing more than 10.0% oxygen.
- B. Not for use in atmospheres containing vapors, gases, or dusts.
- C. Do not exceed maximum use concentrations established by regulatory standards.
- D. Failure to properly use and maintain the product could result in injury or death.
- E. All approved respirators shall be labeled, dated, used, and maintained in accordance with MSDS, OSHA, and other applicable regulations.
- F. Never substitute, modify, add, or use parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- G. Refer to user instructions and/or maintenance manuals for information on use and maintenance of these respirators.
- H. NIOSH does not evaluate respirators for use in high-altitude.



WARNING:

NIOSH APPROVED: N95
95% efficient against particulate aerosols free of oil. This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting disease or infection. Misuse may result in sickness or death.



1. Preshretch top and bottom straps before placing respirator on the face.



2. HOLD the respirator in your hand, with the nosepiece at your fingertips, allowing the headbands to hang freely below your head.



3. Position the respirator under your chin with the nosepiece up. Pull the top strap over your head resting it high at the top back of your head. Pull the bottom strap over your head and position it around the neck below the ears.



4. Place your fingertips from both hands at the top of the metal nosepiece. Using two hands mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece. Ploching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.



5. Perform a User Seal. Check prior to each wearing. To check fit, place both hands completely over the respirator. If air leaks around nose, headband or the straps back along the sides of your head, if you CANNOT achieve a proper fit, DO NOT enter the contaminated area. See your supervisor.

Removal Instructions:

See step 1 of Fitting Instructions and cup the respirator by hand to maintain position on face. Pull bottom strap overhead. Still holding respirator in position, pull top strap over your head and remove respirator.
This respirator contains no components made from natural rubber latex.

NIOSH Approved: N95

All have 95% filtration efficiency against solid and liquid aerosols that do not contain oil.
PROTECTION:
N95-Particulate Filter (95% filter efficiency level) Effective against particulate aerosols free of oil; time use restrictions may apply.

CAUTIONS AND LIMITATION:

- A. Not for use in atmospheres containing less than 19.5% oxygen.
- B. Not for use in atmospheres containing vapors, gases, or dusts.
- C. Make use PPE (OSHA 29 CFR 1910.134) and follow the manufacturer's instructions for use.
- D. Failure to properly use and maintain the product could result in injury or death.
- E. All approved respirators shall be labeled, dated, used, and maintained in accordance with MSDS, OSHA, and other applicable regulations.
- F. Never substitute, modify, add, or use parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- G. Refer to User Instructions and/or maintenance manuals for information on use and maintenance of these respirators.
- H. NIOSH does not evaluate respirators for use in high-altitude.

CH-07-01-A